

# 2014 BioNanotechnology Summer Institute

Cancer Nanotechnology and Cell Mechanics

an NSF-IGERT and  
NIH-CNTC collaboration

July 28–August 8, 2014

University of Illinois at Urbana-Champaign



## MICROFLUIDICS AND ENABLING TECHNOLOGY

### LAB MODULE 1: Fabrication of PDMS-based Microfluidics

Location: 3300 Micro and Nanotechnology Laboratory (MNTL)


Instructors: Angana Senpan, Micro and Nanotechnology Laboratory, BioLab Manager, and Jiwook Shim, Electrical and Computer Engineering

#### Purpose and Expected Outcome:

The purpose of this laboratory module is to provide an introduction and hands-on demonstrations of both the micro-fabrication of PDMS devices and the methods of controlling fluid flow within the device. We will start with a SU-8 master and fabricate devices in PDMS, assemble the chip, and modulate flow while monitoring polystyrene beads and/or food dyes within the device.

#### Overview of Polydimethyl Siloxane (PDMS) Device Fabrication:

PDMS device fabrication is one of the easiest methods for the rapid prototyping of microfluidic devices. The main steps in the fabrication process are sketched in Scheme 1 below.



**Scheme 1.** Schematic overview of PDMS Device Fabrication.

In this experimental module, we will perform the steps outlined above in yellow. The general principles of microfabrication are important considerations for the production of ideal masters. This includes the design and manufacture of positive resists and matching resists to the proper reticles or transparencies. The bonding of the elastomer device will be discussed, a process that is accomplished with an UV-ozone source. Post-bonding surface modifications of PDMS are additional surface chemistry alterations that



## Equipment, Materials, and Supplies:

- Lab coats, gloves, safety glasses
- SU-8 Silicon wafer masters
- PDMS Resin - Dow Corning Sylgard 184 Part A, Part B
- Scale
- Weigh boats
- Stirring bars
- Vacuum jar (desiccator) with vacuum pump (or house vacuum)
- Surgical knives with blades
- Cutting surface (Petri dish lid)
- Sharpened blunt needles to punch holes (See chart at end of document)
- Plasma (oxygen, water, air) source for non-reversible bonding

## Module Outline and Workflow:

In this lab module, participants will get hands-on experience casting (pouring) PDMS over a silicon master device. Because PDMS takes > 2 hours to cure, another set of devices will be prepared ahead of time for cutting.

## Protocol

1. PDMS pouring
  - 1.1. Put on a clean pair of gloves, lab coat, and eye glasses and a face mask.
  - 1.2. Place the clean master mold in a Petri dish. The master can be cleaned of dust or debris that may have accumulated by blowing it with the nitrogen gun.
  - 1.3. On the scale, weigh out and mix PDMS (1:10 ratio) into a weigh boat. Do this by first weighing out 15 g of polymer base, and then add 1.5 g of curing agent, for 16.5 g.
  - 1.4. Please dispose of any extra pre-cured PDMS into a 50 mL conical tube for reuse. PDMS can be stored for overnight at 4°C or for weeks at -20°C without noticeable loss of performance.
  - 1.5. Mix the pre-cured PDMS with a stir bar. Be sure to both swirl and fold the mixture to ensure that the curing agent is evenly distributed.
  - 1.6. Pour the PDMS into the SU-8 master mold placed in a Petri dish.
  - 1.7. Degas the PDMS by placing the mixed pre-cured PDMS in the vacuum desiccator and evacuating the chamber. Bubbles will appear, rise to the surface of the mixture and pop. Degas the mixture for a minimum of 2.0 min. This step may be repeated to completely remove bubbles. Degassing is complete when there are no longer bubbles visible in the mixture. Once all bubbles have been removed, cover the Petri dish and place in an oven at 65-80 °C for 2-6 hrs to cure the PDMS.
2. PDMS release
  - 2.1. Remove the PDMS casting from the oven and place on a clean bench top.
  - 2.2. Using an X-acto knife with a new blade, make straight cuts about 0.5-1 cm from the edge of the master mold. To make each cut, sink the point of the knife vertically into the PDMS until it reaches the polystyrene Petri dish. Keep the knife perpendicular to the master and follow the outline of the master. Make sure to maintain pressure on the knife such that the tip is always in contact with the plastic dish substrate. Continue cutting until the PDMS-master device can be released from the Petri dish with forceps.
  - 2.3. Once all the edges have been liberated, lift the mold up and out of the Petri dish as demonstrated by the instructor. Then carefully peel away the remaining portions of the cured PDMS from the underside (the side without resist features) of the master. Discard this and any excess PDMS.
  - 2.4. Place the released PDMS layer in the lid of your Petri dish with the channel features up.
  - 2.5. With a straight-edge razor, block off the edges of the PDMS to produce a flat PDMS structure that can be bonded to the microscope slide as demonstrated by the instructor.

3. Fluidic port punching
  - 3.1. Align a blunt needle or clean dermal biopsy punch with the first port you will punch.
  - 3.2. Adjust the needle or puncher so that it is as vertical as possible. Push the puncher through the PDMS until you touch the plastic Petri dish. Remove the puncher.
  - 3.3. Push the yellow stick into punched hole and the puncher to drive out the cored section of PDMS.
  - 3.4. Retrieve and discard the cored section from the under side of the device using a pair of forceps. Repeat steps 3.1 to 3.4 for each port.
  - 3.5. Place the punched PDMS device onto a Petri dish with feature side UP.

Once the devices have been poured, cut, and punched, you can *reversibly* bond the PDMS to microscope slides through conformal contact, or irreversibly through plasma treatment or UV-ozone. Equipment varies between labs for cleaning and activating PDMS replicas and glass slides for covalent bonding. Recommended conditions for cleaned PDMS in an oxygen plasma chamber or alternatively with a UV-ozone system that includes (100 mW, 2% oxygen, 35 s) in a PX-250 plasma chamber (March Instruments, Concord, MA). After plasma or UV-ozone treatment, immediately place the oxidized PDMS in contact with the glass to irreversibly bond the surfaces. Chambers are then baked at 70 °C for 10-30 min following bonding.

4. Device bonding (with UV-ozone or oxygen plasma surface treatment)
  - 4.1. Follow protocols specified by the local equipment owner for the oxygen plasma treatment.
  - 4.2. Using forceps or tweezers, place the PDMS device with the feature side facing upwards to be exposed to UV-ozone or oxygen plasma.
  - 4.3. Using forceps or tweezers, place clean glass slides next to the device to be bonded.
  - 4.4. If there is any visible dust particles on the PDMS or the device or slide to be bonded to the PDMS, wipe with a clean lint-free kim wipe or clean room cloth soaked in IPA.
  - 4.5. Place cover on the UV-ozone or oxygen plasma source. For UV-ozone, ensure that the device is approximately 3-5 mm from the UV lamp, which is housed in the cover. Distance requirements are not necessary for oxygen plasma treatments.
  - 4.6. Expose device to oxygen plasma or UV-ozone for 3-5 minutes.
  - 4.7. Remove cover and retrieve PDMS channels using forceps or tweezers, grasp PDMS slab from its side and flip device over onto the glass side so that the features are bonded against the glass.
  - 4.8. Place the devices on a hotplate at 70 °C for 5-10 minutes.

The reactive silanol bonds at the surface of the PDMS will slowly diffuse back into the bulk of the PDMS elastomer. For longer PDMS surface activation and quicker PDMS-substrate bonding, solvent extracted PDMS<sup>5</sup> can be used. Therefore, chemical modification of the PDMS surface should immediately follow the oxygen plasma/ozone bonding for optimal results.

5. Actuating flow control in microfluidic devices can be achieved through a number of methods, during this lab course you will have the opportunity to observe a few of these approaches. They include passive pumping, gravity flow, and syringe pumping.
  - 5.1. Demonstration of passive pumping
  - 5.2. Demonstration of gravity flow
  - 5.3. Demonstration of syringe pumping
  - 5.4. Student group discussions: discuss the advantages of each method of flow control.

## References:

1. McDonald, J. C.; Duffy, D. C.; Anderson, J. R.; Chiu, D. T.; Wu, H.; Schueller, O. J.; Whitesides, G. M., Fabrication of microfluidic systems in poly(dimethylsiloxane). *Electrophoresis* **2000**, 21, (1), 27-40.
2. Makamba, H.; Kim, J. H.; Lim, K.; Park, N.; Hahn, J. H., Surface modification of poly(dimethylsiloxane) microchannels. *Electrophoresis* **2003**, 24, (21), 3607-19.
3. *Silicon Compounds: Silanes and Silicones*. Gelest, Inc.: Morrisville, PA, 2004; p 560.
4. Hermanson, G. T.; Mallia, A. K.; Smith, P. K., *Immobilized Affinity Ligand Techniques*. Academic Press: San Diego, CA, 1992; p 454.
5. Millet LJ, Stewart ME, Sweedler JV, Nuzzo RG, Gillette MU., Microfluidic devices for culturing primary mammalian neurons at low densities. *Lab Chip*. **2007** Aug, 7, (8):987-94.

## Dermal biopsy punches:

Integra\* Miltex\* Standard Biopsy Punches

Assorted sizes, 2mm, 3mm, 4mm, 5mm, and 6mm; 10 each

Fisher Scientific Cat#: 12-460-414

MILTEX INC product no.: 33-38

## Sources for tubing and blunt end needles:

Small Parts, Inc.

13980 N.W. 58th Court

P.O. Box 4650

Miami Lakes, FL 33014-0650

<http://www.smallparts.com>

## Recommended parts descriptions:

Part Description	Usage	Inner Diameter	Outer Diameter	Small Parts Part #
20G x 1/2" Stainless Steel Blunt Needles	Needles for cutting holes	0.023"	0.036"	NE-201PL-C
22G x 1/2" Stainless Steel Blunt Needles	Needles for direct injecting	0.016"	0.028"	NE-221PL-C
30G x 1/2" Stainless Steel Blunt Needles	Needles for tubing	0.006"	0.012"	NE-301PL-C
Tygon Tubing	Connect needles to device	0.01"	0.03"	TGY-010-5C

PDMS – Dow Corning Sylgard 184

<http://www.dowcorning.com/applications/search/products/Details.aspx?prod=01064291&type=PROD>

Available from Ellsworth Adhesives: Part#: 182 SIL ELAST KIT .5KG

UV-Ozone source – PSD-UV, Novascan Technologies



# 2014 BioNanotechnology Summer Institute

Cancer Nanotechnology and Cell Mechanics

an NSF-IGERT and  
NIH-CNTC collaboration

July 28–August 8, 2014

University of Illinois at Urbana-Champaign



## MICROFLUIDICS AND ENABLING TECHNOLOGY

### LAB MODULE 2: Dielectrophoresis

Location: 3300 Micro and Nanotechnology Laboratory (MNTL)

Instructors: Jiwook Shim, Electrical and Computer Engineering, and Angana Senpan, Micro and Nanotechnology Laboratory, BioLab Manager

#### Purpose and Expected Outcome:

The purpose of this laboratory module is to provide an introduction and a hands-on demonstration of microfluidic dielectrophoresis (DEP). The DEP devices are electrodes patterned on a printed circuit board (PCB) that are brought into contact with a very thin glass coverslip that is attached to PDMS microfluidic channels. This assembly will be used to demonstrate trapping and concentration of micro-particles. The students will be able to vary DEP waveform characteristics and related experimental parameters to examine the interactions between the particles and the electric fields generated by the inter-digitated electrodes patterned on the chip. The expected outcome is for the students to gain a basic understanding of dielectrophoresis and its potential applications for biology and medicine.

#### Overview of Dielectrophoresis:

Dielectrophoresis is the electrokinetic movement of electrically polarizable particles in non-uniform electric fields. The non-uniform electric field exerts a force to each end of the polarized particle, with the difference in the magnitude of the two forces determining the direction of particle mobility. DEP occurs for charge-neutral particles and for both DC and AC excitation. Forces in the direction of increasing electric field strength (positive DEP) occur when the permittivity of the particle ( $\epsilon_p$ ) exceeds that of the medium ( $\epsilon_m$ ), whereas particles are pushed towards lower levels of electric field strength when  $\epsilon_p < \epsilon_m$  (Fig. 1), making the polarity of the applied field irrelevant.




Fig. 1: (a) Schematic illustration of dielectrophoresis and (b) representative p-DEP on a cell.

The DEP force on a homogeneous and lossless ( $\sigma/\omega\epsilon \ll 1$ ) sphere of radius ( $r$ ) is:


$$F_{\text{DEP}} = 2\pi r^3 \epsilon_m \operatorname{Re} \left\{ \frac{\epsilon_p^* - \epsilon_m^*}{\epsilon_p^* + 2\epsilon_m^*} \right\} \nabla |E|^2,$$

where complex permittivity ( $\epsilon^*$ ) is equal to  $\epsilon + \sigma / j\omega$  and  $\operatorname{Re}\{\dots\}$  is the *Clausius-Mossotti* factor. Note that this force scales with  $V^2$  and  $r^3$ . Viscous drag on a spherical particle is described by Stokes flow as  $F_d = 6\pi r \eta u$ , where  $\eta$  and  $u$  are the fluid viscosity and velocity, respectively.


Many biologically important particles are polarizable. For example, cells can be described by a shell model, in which the particle is assumed to be composed of a thin membrane surrounding the core, with specified conductivity and permittivity, allowing the DEP force to be estimated.

The implementation of DEP requires patterning of conductive electrodes for the application of non-uniform electric fields, and preferably a layer of insulator on the electrodes to prevent electro-thermal induced reactions at the electrode interface. Alternatively, the electrodes can be integrated within a microfluidic channel for improved functionality depending on the application of interest. In this experiment, we take advantage of a simple and effective microfabrication process that eliminates direct exposure of the target particles to the electrodes (Fig. 2(a)).

The PCB is widely used in the electronics industry to provide mechanical support for electrical connections within the electronic devices. Due to broad use of PCBs, they can be made-to-order from PCB manufacturers. Furthermore, a biochip composed of a PDMS microfluidic channel and a microscope coverslip (Fig. 2(b)) isolates the sample from the electrodes and reduces the risk of cross-contamination between experiments performed on the same PCB.



**Fig. 2.** (a) A schematic diagram of the conventional approach in implementing DEP. Electrodes and substrate should be disposed of after each experiment. (b) A schematic diagram of the PCB-based DEP implementation. DEP electrodes are easily fabricated on a PCB and can be reused again and again since the electrodes are not in direct contact with the sample.



**Fig. 3:** Numerical analysis of the electric potential generated by the PCB electrodes. Electrodes in red and in black are set to +1V and -1V respectively. *Inset:* Square of the electric field magnitude,  $|E|^2$ , on the dotted line 10μm above the glass cover layer. Black and red bar shows the position of the electrodes. The electric field intensity is weakest at the center of the electrode and strongest between the electrodes.

### Equipment, Materials, and Supplies:

- Printed Circuit Board (PCB) with electrode patterns fabricated on the surface
- Polydimethylsiloxane (PDMS) microfluidic channels and *No.0* thickness (~100 μm) coverslips
- Oxygen plasma system
- Function generator, oscilloscope, BNC T-junction and connectors
- AC voltage amplifier
- Upright Optical Fluorescence Microscope with 10-20X objective, with CCD camera.
- Cells or Poly-styrene beads
- Mineral oil
- Double-sided tape
- Soldering iron and solder
- Dermal tissue biopsy hole puncher.


### Module Outline and Workflow:

The students are expected to gain an understanding of dielectrophoretic (DEP) phenomena during an experiment that involves device assembly, microfluidic connection and variation of the applied voltage waveform necessary for DEP manipulation of beads and/or cells. The following are steps the student will take to investigate DEP.


1. Enable connection of PCB electrodes to BNC connectors by soldering two wires onto ends of PCB board. One wire is for positive polarity, the other is for negative polarity.
2. Create PDMS/glass microfluidics chip
  - a. Punch inlet/outlet ports on both ends of channel with PDMS hole puncher by gentle pressure and twisting (punch into channel side first)
  - b. If channel side of PDMS is noticeably dirty, spray off debris with ethanol and dry with nitrogen (PDMS needs to be clean for glass bonding) or remove any debris using adhesive tape
  - c. Lightly press PDMS (channel side down) onto the coverslip to create a bond.
  - d. Align the pipette's tip to the punched inlet and inject  $\sim 20 \mu\text{L}$  of bead suspension into one of the chip's ports to load the channel.
3. Perform DEP on the injected beads/cells and optimize the input voltages and frequencies.
  - a. Mount PCB board and microfluidic chip assembly on microscope stage with double-sided tape and connect the PCB to the amplifier output using mini-grabber BNC connectors (with signal generator and AC voltage amplifier turned off).
  - b. Vary the input voltage and frequency (1 to 12 MHz) to optimize the focusing ability of the device.

**NOTE:** do not allow the AC voltage amplifier's power output to go above 200 Watts! This may cause the PCB to melt or to damage the amplifier's internal components. Also, do not go below 1 MHz to ensure the aforementioned problem does not happen.

  - c. Measure the time it takes for the particles to align to the pattern during DEP. From this data, estimate particle velocities and calculate the hydrodynamic drag force ( $F_{drag} = 6\pi\eta Rv$ , where  $\eta$  = fluid viscosity = 1mPa's for DI,  $R$  = radius of the particle = 7.5  $\mu\text{m}$  (for polystyrene beads), and  $v$  = velocity of the particle motion), which is equal to the DEP force exerted on the particles. This will help in gaining an intuitive understanding of dielectrophoretic manipulation of particles within microfluidic devices (Fig. 4 and 5).



**Fig. 4:** The device is made from a PCB and a thin glass coverslip. The glass coverslip serves as an insulator between the solution and the electrode pattern.



**Fig. 5.** Dark field images of the fluorescently labeled 3 $\mu$ m polystyrene beads during DEP. The trajectories of the polystyrene beads are visible as short lines, indicating the traveled distance during exposure time of the camera. After applying the AC voltages for 36 seconds, most of the beads become focused on top of the electrodes by negative DEP force

#### **Related References:**

1. P. R. C. Gascoyne, J. Vykoukal, "Particle separation by dielectrophoresis," *Electrophoresis*, vol. 23, pp. 1973-83, 2002.
2. M. P. Hughes, *Nanoelectromechanics in Engineering and Biology*. CRC Press, 2002.
3. H. Li and R. Bashir, "On the design and optimization of micro-fluidic dielectrophoretic devices: a dynamic simulation study," *Biomedical Microdevices*, vol. 6, no. 4, pp. 289-95, 2004.
4. H. Li, Y. Zheng, D. Akin, and R. Bashir, "Characterization and modeling of a micro-fluidic dielectrophoresis filter for biological species," *IEEE/ASME Journal of Microelectromechanical Systems*, vol. 14, no. 1, pp. 105-111, 2005.
5. B. M. Taff and J. A. Voldman, "A scalable addressable positive-dielectrophoretic cell-sorting array," *Analytical Chemistry*, vol. 77, pp. 7976-83, 2005.

